

areaDetector Update: EPICS V4 and More

Mark Rivers

GeoSoilEnviroCARS, Advanced Photon Source

University of Chicago



Outline

- Last TWG areaDetector talk was April 2016
 - Releases since then: R2-5, R2-6 (very soon)
- EPICS V4 driver and plugin
- Other new and improved drivers
- Other new and improved plugins
- ImageJ plugins
- Future plans

EPICS Version 4 in a Nutshell

- New Protocol, “pvAccess”
- Structured data
- Introspection interface, “pvData”
- Dynamic typing
- Streaming
- Standard Scientific Types
- RPC and put/get
- New smart database
- Codec based transport
- All APIs in C++ and Java
- Python and Matlab
- High Performance
- High Reliability

```
$ eget -s XCOR:LI24:900:TWISS
non-normative type
structure
  double energy 5.00512
  double psix 37.7625
  double alphax 13.6562
  double betax -2.78671
  double etax -0.00698294
  double etaxp 0.00107115
  double psiy 31.9488
  double alphay 116.762
  double betay 5.2592
  double etay 0
  double etayp 0
```

Figure: Example EPICS v4 pvAccess method “eget”, getting PV whose value is a structure giving the Courant-Snyder parameters at a corrector magnet. This example is from the LCLS control system at SLAC.

The EPICS V4 “Normative Types”

The Normative Types Spec defines a standard for commonly used data types, <http://epics-pvdata.sourceforge.net/alpha/normativeTypes/normativeTypes.html>

5. General Normative Types

1. NTScalar
2. NTScalarArray
3. NTEnum
4. NTMatrix 
5. NTURI 
6. NTNameValue
7. NTTTable 
8. NTAttribute

```
$ eget -s XCOR:LI24:900:RMAT
  0.0727485    0.0289316          0          0    0.0652488    0.00125391
  0.0578214    0.0391775          0          0    -0.027185   -0.000192344
                0          0    0.00943029    1.14291          0          0
                0          0   -0.0013367   -0.0348832          0          0
-0.000370971 -0.000283933          0          0   -0.0182387  -0.000198345
  0.10031     0.018722          0          0   -10.5721    -0.179568
```

```
$ eget pva://mccas0.slac.stanford.edu:39633/QUAD:LTU1:880:RMAT?type=design
```

```
$ eget -s LCLS:ELEMENTS
ELEMENT      ELEMENT_TYPE      EPICS_DEVICE_NAME      S_DISPLAY      OBSTRUCTION
CATHODE      MAD                CATH:IN20:111          2014.7         N
SOL1BK      MAD                SOLN:IN20:111          2014.7         N
CQ01        MAD                QUAD:IN20:121          2014.9         N
SOL1        MAD                SOLN:IN20:121          2014.9         N
XC00        MAD                XCOR:IN20:121          2014.9         N
...
```

6. Specific Normative Types

1. NTMultiChannel
2. NTNDArray 
3. NTContinuum
4. NTHistogram
5. NTAggregate





Figure: An extract of the Table of Contents of the Normative Types Specification document, together with examples of 4 selected types


EPICS Version 3 basic block diagram

EPICS in the nominal usage: An EPICS client communicates over Channel Access (CA) protocol to an Input/Output Controller (IOC)
Channel Access server (module rSrv in an IOC)



EPICS Version 4 is an extension of V3

V4 IOC == V3 IOC + pvAccess Server




Use Case: Network efficient acquisition of archived meta data

Presently, only 1 PV per pvAccess channel. But plan is to get/monitor a group of PVs through one pvAccess channel.


EPICS Version 4 includes CA

The pvAccess API includes Channel Access support, **so one client lib does both**




EPICS Version 4 new database

A new smart database, “pvDatabase” can be used for data assembly and processing




EPICS Version 4 middleware support

RPC and Service Oriented Architecture (SOA)




NDPluginPva

- New plugin that converts NDArrays into the EPICSv4 normative type NTNDArray
- An embedded EPICSv4 server serves the new NTNDArray structure as an EPICSv4 PV
- High performance, ~3.2GB/s shown here
- Can be received by any EPICSv4 client
 - Java, Python, C++ versions of pvAccess
 - CSS has a widget that can display NTNDArrays
 - caQtDM has a new camera widget to display NTNDArrays
 - Tim Madden is writing a version of the Java areaDetector ImageJ plugin based on pvAccess; will talk to this plugin
 - Can include an NTNDArray receiver in another IOC



pvAccess Driver

- New pvAccess driver receives NTNDArrays over the network, converts to NDArrays and calls plugins
- Can be used to run areaDetector IOC and plugins on another machine or in another process
- High performance:
 - ~1.2 GB/s shown here with interprocess communication
 - Saturating 10 Gb Ethernet links has been demonstrated



pvAccess: Command Line Tools

```
>pvinfo 13SIM1:Pva1:Image
CHANNEL : 13SIM1:Pva1:Image
STATE   : CONNECTED
ADDRESS : 164.54.160.82:5075
epics:nt/NTNDArray:1.0
  union value
    boolean[] booleanValue
    byte[] byteValue
    short[] shortValue
    int[] intValue
  ...
  uint[] uintValue
  ulong[] ulongValue
  float[] floatValue
  double[] doubleValue
  ...
  int uniqueId
  time_t dataTimeStamp
    long secondsPastEpoch
    int nanoseconds
    int userTag
  ...
```

```
>pvget 13SIM1:Pva1:Image -r
'field(uniqueId,dimension)'
13SIM1:Pva1:Image
structure
  int uniqueId 886580
  dimension_t[] dimension
    dimension_t
      int size 1024
      int offset 0
      int fullSize 1024
      int binning 1
      boolean reverse false
  dimension_t
    int size 1024
    int offset 0
    int fullSize 1024
    int binning 1
    boolean reverse false
```

NDPluginPVA and pvAccess Driver Demo

NDDriverStdArrays

- Driver that allows any EPICS Channel Access client to create NDArrays in an areaDetector IOC.
- Logical inverse of NDPluginStdArrays
 - That plugin converts NDArrays in an IOC into standard EPICS waveform records for use by Channel Access clients.
 - Also writes to additional records to describe the array structure. (SizeX, SizeY, ColorMode, etc.)
- NDDriverStdArrays
 - This driver receives EPICS waveform records from CA clients
 - Converts them to an NDArray in the IOC.
 - Clients must also write to additional EPICS records to define the structure of the NDArray.
- Similar in concept to pvaDriver
 - pvaDriver only works with EPICS V4 PVAccess clients
 - NDDriverStdArrays works with EPICS V3 Channel Access clients.

NDDriverStdArrays Use Case

David Vine, LBNL

- Display data from sscan record
- Simple Python client reads data from sscan record, writes to NDDriverStdArrays
- Allows real-time display of 2-D scans
- Can write data in standard HDF5 file format
- Any other AD plugins can be used on the scan data

NDDriverStdArrays

- Driver receives 1-D arrays, normally from an EPICS waveform record.
- NDimensions record controls actual number of dimensions to use when creating the NDArrary.
- Dimensions array record controls actual size of each of the NDimensions dimensions.
- ColorMode record can be used to define the color mode (Mono, Bayer, RGB1, RGB2, or RGB3) of the NDArrary.
- Data type of the NDArrary set by the DataType record.
 - Int8, UInt8, Int16, UInt16, Int32, UInt32, Float32, or Float64.
- Datatype (FTVL) of the EPICS waveform record can be
 - CHAR, UCHAR, SHORT, USHORT, LONG, ULONG, FLOAT, or DOUBLE.
- Driver does any required type conversion between the waveform record data and the NDArrary data.

NDDriverStdArrays

- 2 modes of operation, AppendMode=Disable/Enable
 - Disable: each write to the waveform record generates a complete NDArrary.
 - Enable: waveform record contains just part of the NDArrary data; data appended by multiple writes to the waveform record to construct the complete NDArrary. sscan record client would work in this mode
- Records in Enable mode:
 - NewArray, NextElement, Stride, ArrayComplete
 - CallbackMode record controls when the driver does callbacks
 - OnUpdate, OnComplete, OnCommand

Area Detector Control - 13NDSA1:cam1:

Setup

asyn port **NDSA**
 EPICS name **13NDSA1:cam1:**
 Manufacturer **NDDriverStdArrays**
 Model **Software Detector**
 Serial number **No serial number**
 Firmware version **No firmware**
 SDK version **1.1.0**
 Driver version **1.1.0**
 ADCore version **2.6.0**

Connected

Connection

Debugging

Plugins

Readout

Dimensions
 Array size
 # elements
 Next element
 Stride
 Fill value
 Actual size
 Image size (bytes)
 Data type **Float32**
 Color mode **Mono**

Attributes

File

Collect

Images
 # Images complete
 Image mode

Collecting

Acquire

Detector state

Status

Image counter

Image rate

Array callbacks

Buffers

Buffers max/used
 Buffers alloc/free
 Memory max/used (MB)
 Buffer & memory polling

Other

Append mode
 Callback mode

NDDriverStdArrays Demo

Runs Python Unit Tests

- Complete NDArrays
- Line scans (appends line at a time)
- Point scans (appends pixel at a time)
 - 4 detectors signals, sent to 4 ROIs

Pilatus Driver Update

- Now supports Energy and ThresholdEnergy independently
 - ThresholdEnergy is used to set per-pixel discriminators
 - Energy is used to compute the correct flat-field normalization
- Previously EPICS only exposed ThresholdEnergy
 - Camserver was supposed to set $\text{Energy} = 2 * \text{ThresholdEnergy}$
 - Bug in some versions of camserver, e.g. our new 1M CdTe detector. It was freezing Energy at 60 keV, and we collected an entire run of improperly normalized data
- There are occasions when $\text{ThresholdEnergy} = \text{Energy}/2$ is not the optimum setting, i.e. to avoid strong fluorescence.

Pilatus Driver Update

- If camserver is saving TIFF files then driver now reads the TIFFImageDescription tag from the TIFF file.
 - This is a long string that camserver writes to the file containing all of the detector settings, including threshold, energy, etc.
 - Driver adds this information to the NDArray using an NDAttribute called TIFFImageDescription.
- The NDFileTIFF plugin in ADCore R2-6 was changed to write this complete attribute to the TIFFImageDescription tag in the new TIFF file.
- It will also be written by the NDFileNetCDF, NDFileHDF5, and NDFileNexus plugins.
 - However these plugins are currently limited to 256 character string attributes, so some of the information will be lost because the string is longer than 256 characters.

ADAndor

- Added support for Electron Multiplying (EM) Gain.
- Added ability to set the BaselineClamp in the Andor SDK.
- Enforce minimum values of ADShutterOpenDelay and ADShutterCloseDelay based on query of SDK.
- Implemented ReverseX and ReverseY.
- Fixed bug with AndorPreAmpGain; previously it was not actually calling SetPreAmpGain().
- Added support for Full Vertical Binning (FVB) readout mode.
- Added support for EPICS shutter control.

Pixirad driver update

- CdTe pixel array detector available with multiple chips.
 - Detector pool has Pixirad-2, 402x1024. Now working at 1-ID.
- Updated driver to work with new PIII chip.
 - Square 62 um pixels pixels, not hexagonal.
 - Single chip is 402x512 pixels.
 - Pixel mode (PM). Each pixel is independent.
 - Highest rate capability ($> 1\text{MHz}$), no sharp color separation, multiple counts if threshold is low.
 - Neighbor Pixel Inhibit mode (NPI).
 - Only one counter per event allowed to count. The hit is allocated to the pixel receiving the highest fraction of the total charge. Highest position resolution and the lowest possible noise but sharpness of the color separation not completely solved yet
 - Pixel Summing mode (PSM).
 - The signals of 4 neighbor pixels are summed together in each pixel to correctly evaluate the total energy of any event involving up to 4 pixels.

Other detectors

- New ADCameraLink, ADPCO drivers
 - Drivers for PCO cameras from Tim Madden
- New ADLambda driver
 - Driver for Lambda pixel array detectors from John Hammonds
- Many detectors (Pilatus, Prosilica, Point Grey, marCCD, mar345, Andor, etc.)
 - Added support for SerialNumber, FirmwareVersion, SDKVersion, DriverVersion, and ADCoreVersion which were added in ADCore R2-6.
- ADSimDetector, ADCSimDetector, pvaDriver
 - Moved out of ADExample and into their own repositories
- ADExample no longer used

ADSupport

- New repository ADSupport
- Source code for all 3rd party libraries used by ADCore.
 - hdf5, jpeg, netCDF, nexus, szip, tiff, xml2, zlib
- ADBinaries no longer used.
 - Pre-built libraries for Windows too difficult to maintain, now all built from source using EPICS build system
- ADSupport can be built for
 - Windows (Visual Studio or MinGW, 32/64 bit, static/dynamic),
 - Linux (currently Intel architectures only, 32/64 bit),
 - Darwin
 - vxWorks (currently big-endian 32-bit architectures only).
- vxWorks previously supported only netCDF file plugin
 - All file saving plugins now supported on vxWorks 6.x (TIFF, JPEG, netCDF, HDF5, Nexus).
 - HDF5 and Nexus are not supported on vxWorks 5.x because the compiler is too old.
 - Could imagine migrating saveData to HDF5 since it now works on all platforms.

Optional Support

- All 3rd party libraries now optional.
 - For each library XXX 4 Makefile variables control support for that library. XXX can be JPEG, TIFF, NEXUS, NETCDF, HDF5, XML2, SZIP, and ZLIB.
 - WITH_XXX YES then drivers or plugins that use this library will be built. NO then drivers and plugins that use this library will not be built.
 - XXX_EXTERNAL YES then the library is not built in ADSupport, but assumed to be found external to the EPICS build system. NO then the XXX library will be built in ADSupport.
 - XXX_DIR If defined and XXX_EXTERNAL=YES then the build system will search this directory for the XXX library.
 - XXX_INCLUDE If defined then the build system will search this directory for the include files for the XXX library.
- All EPICS modules except base and asyn are now optional.

NDPluginOverlay

- New Ellipse shape to draw elliptical or circular overlays.
- Improved efficiency by only computing the coordinates of the overlay pixels when the overlay definition changes or the image format changes.
 - Important for ellipse because it uses trig functions
- Added CenterX and CenterY parameter for each overlay.
 - Can now specify the overlay location either by PositionX and PositionY, (upper left corner), or by CenterX and CenterY,
 - If CenterX/Y is changed then PositionX/Y will automatically update, and vice-versa.
- Changed the meaning of SizeX and SizeY for the Cross overlay shape.
 - Previously the total size of a Cross overlay was $\text{SizeX} * 2$ and $\text{SizeY} * 2$.
 - It is now SizeX and SizeY.
 - Now consistent with the Rectangle and Overlay shapes, i.e. drawing each of these shapes with the same PositionX and SizeX/Y will result in shapes that overlap in the expected manner.

Other Plugin Enhancements

- All plugins: ElapsedTime PV for last execution.
 - This gives the execution time in ms the last time plugin ran.
 - Works both with BlockingCallbacks=Yes and No.
 - Convenient for measuring plugin performance without having to run detector at high frame rates.
- NDPluginROI
 - Added CollapseDims to optionally collapse (remove) output array dimensions whose value is 1. For example an output array that would normally have dimensions [1, 256, 256] would be [256, 256] if CollapseDims=Enable
- NDPluginStats
 - Added calculations of 3rd and 4th order image moments, this provides skewness and kurtosis.
 - Added eccentricity and orientation calculations.
- NDPluginFile
 - If the NDArray contains an attribute named FilePluginClose and the attribute value is non-zero then the current file will be closed.

Base Class Enhancements

- `asynNDArrayDriver`, `NDArrayBase.template`, `NDPluginBase.adl`, `ADSetup.adl`, all plugin adl files
 - `ADCoreVersion_RBV`, `DriverVersion_RBV`.
 - Show the version of `ADCore` and of the driver or plugin that the IOC was built with.
- `ADDriver`, `ADBase.template`, `ADSetup.adl`
 - `SerialNumber_RBV`, `FirmwareVersion_RBV`, and `SDKVersion_RBV`.
 - Show the serial number and firmware version of the detector, and the version of the vendor SDK library that the IOC was built with.
- Because `NDPluginBase.adl` and `ADSetup.adl` grew larger all plugin and driver adl files have changed their layouts.

Base Class Enhancements

NDPva.adl@corvette.cars.aps.anl.gov

13SIM1:Pva1:

asyn port PVA1
 Plugin type NDPluginPva
 ADCore version 2.6.0
 Plugin version 2.6.0
 Array port SIM1 SIM1
 Array address 0 0
 Enable Enable
 Min. time 0.000 0.000
 Callbacks block No
 Queue size/free 20 30
 Array counter 0 326820
 Array rate 784.00
 Execution time 0.142 msec
 Dropped arrays 0 0
 # dimensions 2
 Array Size 1024 1024 0
 Data type Float32
 Color mode Mono
 Bayer pattern RGGB
 Unique ID 326820
 Time stamp 854908701.126
 Attributes file
 Array callbacks Enable
 asyn record

EPICS V4 PV Name
 V4 PV Name 13SIM1:Pval:Image

Andor.adl@corvette.cars.aps.anl.gov

Andor Detector Control - 13ANDOR1:cam1:

Setup

asyn port ANDOR
 EPICS name 13ANDOR1:cam1:
 Manufacturer Andor
 Model DU401_BU2
 Serial number 13278
 Firmware version 3.255
 SDK version 2.99.30001.0
 Driver version 2.5.0
 ADCore version 2.6.0
 Connected
 Connection
 Debugging

Shutter

Shutter Type
 Andor Shutter Mode
 External shutter
 Status: Det. EPICS
 Open/Close
 Delay: Open 0.010 Close 0.020
 EPICS shutter setup

Collect

Exposure Time 1.000 1.000
 Accumulate Period 2.000 2.453
 Acquire Period 5.000 2.453
 # Accums/Image 1 1
 # Exposures Complete 1
 # Images/Acquis. 10 10
 # Images Complete 1
 Image Mode Single
 Trigger Mode
 Readout Mode Image
 Done
 Acquire
 Detector State Idle
 Detector Status IDLE. Waiting on instructions
 Andor Message
 Image Counter 0 45
 Image Rate 0.00
 Array Callbacks Enable

Plugins

Readout

	X	Y
Sensor Size	1024	127
Binning	1	1
Region Start	0	0
Region Size	1024	127
Reverse	<input type="checkbox"/> No	<input type="checkbox"/> No
Image Size	1024	127
Image Size (bytes)	260096	
Pre-amp Gain	1.00	1.00
EM Gain	0	0
EM Gain Mode	<input type="text" value="8 bit DAC"/> 8 bit DAC	<input type="text" value="8 bit DAC"/> 8 bit DAC
EM Gain Adv.	<input type="text" value="Disabled"/> Disabled	<input type="text" value="Disabled"/> Disabled
ADC Speed	<input type="text" value="0.10 MHz"/> 0.10 MHz	<input type="text" value="0.10 MHz"/> 0.10 MHz
Data Type	<input type="text" value="UInt16"/> UInt16	<input type="text" value="UInt16"/> UInt16

Attributes

File

File

Driver File I/O

Cooler

Cooler On
 Temperature -20.0 -19.9
 Status Stabilized at set point

PAL File

Path

Viewers

- EPICS_AD_Viewer enhancements
 - Automatically set the contrast when a new window is created.
 - Eliminates the need to manually set the contrast when changing image size, data type, and color mode in many cases.
 - When image window is automatically closed and reopened because the size, data type, or color mode changes the new window is now positioned in the same location as the window that was closed.
- Gaussian Profiler
 - Like the normal ImageJ line profile in “live” mode, but with fitting of Gaussian and live display of fit parameters

Viewers

- EPICS_AD_Controller. Allows using the ImageJ ROI tools (rectangle and oval) to graphically define the following:
 - The readout region of the detector/camera
 - The position and size of an ROI (NDPluginROI)
 - The position and size of an overlay (NDPluginOverlay)
 - The plugin chain can include an NDPluginTransform plugin which changes the image orientation and an NDPluginROI plugin that changes the binning, size, and X/Y axes directions. The plugin corrects for these transformations when defining the target object.
 - Chris Roehrig wrote an earlier version of this plugin.

Viewer Demo

Future Ideas: ADCore R3-0

- Put more functionality into ADDriver base class
 - Currently it does not do much, all code is in each driver for:
 - Doing callbacks to plugins
 - Processing new exposure time with writeFloat64 function
 - writeFloat64 in ADDriver base class would call setExposure() in derived class
 - Derived class would call ADDriver::doPluginCallbacks(), which would handle setting attributes, getting timestamp, calling plugins, etc.
- Model 3 motor driver and quadEM, which also use asynPortDriver, are written this way
- Demultiplexor/multiplexor plugin
 - Allow multiple plugins to work on the same data stream when it saturates a single core
- Simplify file saving modes (no more Single, Capture, Stream)

Future Ideas: ADCore R4-0

- Change NDArray to NTNDArray for passing data to plugins
- Utilize EPICS V4 infrastructure
- Smart pointers automatically eliminate all unnecessary copying
- V4 clients can immediately receive data with no need to convert to waveform records
- Bruno Martins has demonstrated this working
- Seems like a good path for the future