

2. Spec commands: local commands at 4idd

General commands:

- qdo macro1.mac
 - load “macro1.mac” located in the current directory
- ascan, dscan
 - absolute scan, relative scan
 - alternatively “lup” instead of “dscan”
- Escan
 - Monochromator energy scan, scan high to low
- pl_xMIN, pl_xMAX, CEN
 - x-positions of minimum, maximum, center of mass from the last plot
- splot_abs, splot_dabs
 - replot absorption ($\ln(-\text{Det}/\text{Mon})$) to screen
 - replot derivative of absorption
- cplot_abs, cplot_dabs: send plots to printer
- cplot_plot: plot last scan to printer

2. Spec commands: local commands at 4idd

General commands: counters

- counters
 - Sets Monitor and Detector
- ct [uct]
 - Counts for specified time
- plotselect
 - Selects multiple detectors to plot during scan

```
6134.FOURC4ID> counters
Current counter configuration:
Num          Name  Mnemonic
0           Seconds  sec
1            IC1   ic1
2            IC2   ic2
4            IC4   ic4    <- Assigned to MON
5            IC5   ic5    <- Assigned to DET
21           Iaps   iaps
22      SplitSum  spsum
23      SplitPos  sppos

Counter number for monitor, -1 disables (4)?
Counter number for detector (5)?

6135.FOURC4ID> █
```

2. Spec commands: local commands at 4idd

General commands: motors

- mv [umv]
 - Moves a motor to a position [with updating]
 - umv th 0.1
- mvr [umvr]
 - Moves a motor by a relative amount
 - umvr th 0.1 (increase th by 0.1 Deg)
- wa
 - Displays all motor positions
- wm
 - Shows position of specific motor
- lm
 - Sets limits for motor

2. Spec commands: local commands at 4idd

General commands: undulator tracking, offset

- und_setup
 - sets parameters for undulator tracking
- und_tracking_on
 - enables undulator tracking as the monochromator energy is changed.
- Und_Off
 - parameter for undulator offset
 - in keV
 - $E(\text{undulator}) - E(\text{mono})$

```
6147.FOURC4ID> und_setup
Setup for undulator: ID04ds:
  Monochromator Undulator energy offset: (0.2)?
  Keep undulator taper (mm) fixed? (YES)?
  Undulator taper (0-5.0mm): (0)?
6148.FOURC4ID> und_tracking_on
Wed Apr 19 19:29:05 2017. Undulator tracking is on.
6149.FOURC4ID> p Und_Off
0.2
6150.FOURC4ID>
```

2. Spec commands: local commands at 4idd

General commands: Energy scan using qxscan

- qxscan
 - scan as a function of incident energy
- example
 - qxscan 11.215 -1e5
 - first variable is absorption edge
 - second is for counting time
 - counting time<0: monitor count
 - counting time>0: seconds
- qxscan_setup
 - set parameters used in qxscan
 - Pre-edge region
 - Post-edge region (“k” value)

```
-----
6133.FOURC4ID> qxscan_setup
Enter number of pre-edge regions (1)?
Defining pre-edge region 1
Enter energy relative to edge (in eV) (-20)?
Enter energy increment (in eV) (1)?
Enter collection time in seconds (1)?
Enter starting edge region energy (in eV relative to edge): (-15)?
Enter final edge region energy (in eV relative to edge): (15)?
Enter energy increment in edge region: (0.5)?
Enter collection time in seconds: (1)?
Enter number of post-edge regions (1)?
Your edge region ends at k = 1.984610
Defining post-edge region 1
Enter k end (3)?
Enter k increment (0.07)?
Enter collection time in seconds (1)?

Emin = 11.195, Emax = 11.250; Total qxscan points = 81
6134.FOURC4ID>
```

E-E ₀	k
1	0.51
5	1.15
10	1.62
15	1.98
20	2.29
25	2.56
30	2.81
50	3.62
100	5.12

$$k = (2m(E-E_0)/h^2)^{1/2}$$

$$k = (0.2625 \times [E-E_0])^{1/2}$$

2. Spec commands: local commands at 4idd

General commands: Changing incident energy

- moveE
 - Moves the incident energy (in keV)
 - Undulator and monochromater
- dmoveE
 - Moves phase plate additionally
- getE
 - Displays current incident energy

2. Spec commands: local commands at 4idd

Dichroism related

- `dichro_setup`:
 - Set up parameters related to phase plate control in dichroism scans
- `pr_track_on`, `pr_track_off`
 - option to adjust phase plate angle as the incident energy is changed.
- `dmoveE`
 - change incident energy, tracking phase plate angle
 - undulator, monochromator, and phase plate changes in sequence.

```
6131.FOURC4ID> dichro_setup
Dichroism looping schemes:
0: Helicity Switching      <- Current Loop
1: Magnet Switching
2: Helicity (PZT) & Magnet Switching

1) Dichroism looping scheme (0)?
2) Use ABBA switching (No -> AB only) (YES)?
3) Plot flipping ratio (No -> difference) (YES)?
4) Autoset phase retarder offset value (NO)?
4) Settling time after dichro_move (0.05)?

Phase retarder setup:
1) Use Phase retarder #1 pr1th (NO)?
2) Use Phase retarder #2 pr2th (YES)?
4) Use motors for switching (No -> pzts) (NO)?
8) Use Phase retarder #3 pr3th (NO)?

Phase Retarder 2:
Offset amount: (-0.012)?
Epics PV for PR2 PZT: (4idb:E665:2:DC_set_microns.VAL)?
PR2 PZT center (where diamond tweaked up) (15)?
PR2 PZT micron to degrees conversion: (0.0019324)?

6132.FOURC4ID> █
```

3. list of commonly used SPEC commands at 4idd

Screen plot macros

- splot_abs** Plots $-\log(\text{DET}/\text{MON})$ of last scan on the screen
- splot_dabs** Plots the derivative of $-\log(\text{DET}/\text{MON})$ of last scan on the screen
- splot_deriv** Plots the derivative of last scan on the screen

Printer plot macros

- cplot_plot** Prints a plot of the last scan
- cplot_abs** Prints absorption plot ($-\log(\text{DET}/\text{MON})$) of last scan
- cplot_deriv** Prints plot of the derivative of last scan
- cplot_dabs** Prints a plot of the derivative of the absorption

Slit macros

- showslits** Displays current slit positions
- setslit** Sets slit size and position
- slitscan** Scans a slit blade, center, or position
- conflits** Defines which slit to control in SPEC

Filter Macros

- filter** Sets filter position
- Showfilters** Displays current filter positions

Temp Control Macros

- te** Shows temperature of set setpoint
- changetemp** Changes temperature to value (waits before next command)
- tempscan** Scans temperature
- tempscan_6T** Scans temperature with 6 Tesla magnet lakeshore
- ramp** Turns on ramp or Sets ramp rate
- heater** Sets heater range [0/1/2/3/4/5]
- save_temp_on[off]** Saves temperature in data file
- lakeshore1[2/3]** Selects lakeshore controller

Undulator Macros

- moveundE** Moves undulator energy (with backlash)
- undscan** Scans the undulator in energy
- undscangap** Scans undulator gap
- und_setup** Gives offset for energy moves
- und_tracking_on[off]** Turn on[off] undulator tracking of mono

3. list of commonly used SPEC commands at 4idd

Analyzer Macros

anal_setup	Defines analyzer d spacing and type
anal_on[off]	Sets analyzer track in energy
asetE	Sets analyzer ath and atth positions for given energy
calc_anal	Shows ath and atth positions for given energy
go_sigma[pi]	Moves polarization type analyzer sigma or pi position

Disable Macros

disable[enable]	Disable/enable a motor of counter. Use ? to show
------------------------	--

Dichroism Macros

dichro_setup	Sets dichro scan parameters
pr_track_on(off)	Phase retarder track in E (on or off)
dsetEpr1	Sets phaseR1 theta for a given energy
dsetEpr2	Sets phaseR2 theta for a given energy
dmoveE	Moves monochromator and phaseR to a given E
adichro	1 motor dichro scan at a fixed energy
a2dichro	2 motor dichro scan at a fixed energy
Edichro	Energy dichro scan
Edichro2	Energy dichro scan with 3 energy intervals
qxdichro	Energy dichro scan with multiple energy intervals
kepdichro	Dichro scan with magnet current using kepko powersupply
amidichro	Dichro scan with magnet current using ami powersupply
tempdichro	Dichro temperature scan at a fixed energy
cplot_dichro_flip	Prints a plot of signal and flipping ratio
cplot_dichro_diff	Prints a plot of signal and difference
lock (scan)	Starts a lockin scan
cplot_lock	Prints a plot of the DC and lockin signals
pzt_osc_on(off)	Turns on (off) PZT oscillations

Kepko Macros (Electromagnets)

kepko_on(off)	Turns on(off) kepko powersupply
kepkoI	Sets kepko current
kepkoV	Sets kepko voltage
kepscan	Scans kepko current

3. list of commonly used SPEC commands at 4idd

Ami Macros (4Telsa, 6Telsa superconducting magnets)

amifield Sets magnetic field (in kG)
amiscan Scans fields using ami powersupply

EXAFS Macros

qxscan Energy scan with fixed q steps above the edge
qxscan_setup Defines qxscan parameters, related to the edge energy

Reciprocal space Macros

wh Shows current reciprocal space position
ca Calculates motor positions for H K L
br, ubr Moves to H K L reciprocal space position
or0, or1 Defines primary secondary HKL reflections
setlat Sets lattice parameters
setmode Sets conversion mode
uan Moves theta and two theta to position
th2th theta/2theta relative scan
hscan,kscan,lscan, hklscan Reciprocal space scan